INSTRUCTIONS FOR USE

# SR Range

## Digital Handheld Doppler

Instructions for Use · Mode d'emploi · Gebrauchsanleitung · Istruzioni per l'uso ·  
Instrucciones de uso · Instruções de utilização · Bruksanvisning · Gebruiksaanwijzing ·  
Brugervejledning · Bruksanvisning · Käyttöohjeet · Návod k použití

# HUNTLEIGH

A MEMBER OF THE ARJO FAMILY

# Contents

1. Introduction
  - 1.1 Unpacking / Preliminary Checks
2. Safety
  - 2.1 Warnings & Cautions
  - 2.2 Patient Applied Parts
  - 2.3 Intended Use & Indications
  - 2.4 Contraindications
  - 2.5 Patient Population
  - 2.6 Clinical Benefit
3. Warranty & Service
  - 3.1 Service Life
  - 3.2 Maintenance and Repair
4. Product Identification
  - 4.1 Product Controls
  - 4.2 Symbol Identification
  - 4.3 Display Status Bar
5. Prepare the Doppler for Use
  - 5.1 Battery Insertion / Replacement
  - 5.2 Probe Connection
  - 5.3 Change the Doppler Settings
6. Operation
  - 6.1 SRX Probe options
  - 6.2 SR2 & SR3 Waterproof Probes
  - 6.3 Obstetric Monitoring
  - 6.4 Frozen Trace
  - 6.5 Battery Charging
  - 6.6 Data transfer to a PC
  - 6.7 After Use
7. Care and Cleaning
  - 7.1 General Care
  - 7.2 Cleaning and Disinfecting the Doppler
  - 7.3 Cleaning and Disinfecting Probes
8. Troubleshooting
9. Specifications
  - 9.1 Equipment Classification
  - 9.2 FHR Performance
  - 9.3 General
  - 9.4 Environmental
  - 9.5 Standards Compliance
10. Electromagnetic Compatibility

# 1. Introduction

The SR range are multi-function, battery powered, hand-held Dopplers intended for obstetric use. The SRX is compatible with Huntleigh's interchangeable probes. The range provides an audible fetal heart sound and a numeric/graphical display of fetal heart rate. This equipment is for use only by suitably qualified healthcare practitioners and is not intended for use by the patient.

Before using this equipment, study this manual carefully and familiarise yourself with the controls, display features and operation.

Experience with use of ultrasonic dopplers is preferable, but for novice users training material is provided with the online documents. Exposure to ultrasound should be kept As Low As Reasonably Achievable - (ALARA guidelines).

Scan the QR code on the rear cover of this IFU with a smartphone, or visit the Huntleigh website for electronic copies of user literature. All documents are available to download as PDF files. To read them, you must have a PDF reader installed on your device. Alternatively paper copies are available upon request.

## 1.1 Unpacking / Preliminary Checks

On receipt of your Doppler, check that all items are present and undamaged. If items are missing or have been damaged in transit, inform Huntleigh Healthcare immediately.

Digital Handheld Doppler	IFU (this document)	Batteries
Charger & USB lead*	Ultrasound Gel	Carry bag

\*selected Dopplers

# 2. Safety

## 2.1 Warnings & Cautions

- Dopplers are screening tools to aid the healthcare professional. If there is any doubt as to fetal status or suspected errors in reading fetal heart trace, further investigations should be undertaken immediately using alternative techniques.
- Always ensure that Maternal HR or any Artefacts are not misinterpreted as fetal heartrate. If in doubt, feel the mother's pulse during the examination.
- Ensure adequate cleaning of devices to minimise the risk of cross contamination and infections.
- Precautions should be taken when using ultrasound in the presence of pacemakers and defibrillators.
- These devices are not suitable in utero or in conjunction with intravenous contrast.
- Do not use in the presence of flammable gases.
- Do not use in a sterile field unless additional barrier precautions are taken.
- Do not sterilise the product or its accessories. The product will be damaged.
- Do not expose to excessive heat, including prolonged exposure to sunlight.
- Do not dispose of batteries in fire as this can cause them to explode.
- The Doppler is not waterproof and must not be immersed.
- This product contains sensitive electronics, which are susceptible to interference, this will be indicated by unusual sounds.
- Any equipment connected to the USB port must be compliant with IEC 60601-1.
- This equipment must not be modified.

## 2.2 Patient Applied Parts

As defined in IEC 60601-1, the patient applied parts of the SR Doppler are the ultrasound probes.

## 2.3 Intended Use & Indications

The Dopplers are intended for use by qualified healthcare practitioners in primary, acute and community healthcare environments, for the assessment of fetal heart rate.

They are indicated for routine screening of pregnant women of all ages from early gestation through to full term, and for low risk labour management.

## 2.4 Contraindications

- Do not use on broken or fragile skin.
- Do not use on the eye.

## 2.5 Patient Population

The SR Range is suitable for use on all patient populations.

## 2.6 Clinical Benefit

The clinical benefits include but are not limited to:

- screening of the fetal heart rate at all stages of pregnancy from as early as 8-10 weeks, right through to labour
- enabling real-time assessment of the rate and rhythm of the prenatal heartbeat to ensure the baby is developing as expected during the pregnancy.

## 3. Warranty & Service

Huntleigh Healthcare Diagnostic Products Division standard terms and conditions apply to all sales. A copy is available on request. These contain full details of warranty terms and do not limit the statutory rights of the consumer.

**Service Returns:** To return the Doppler, please:

- Clean the product following the instructions in this manual.
- Pack it in suitable packing.
- Attach a decontamination certificate (or other statement declaring that the product has been cleaned) to the outside of the package. (Huntleigh Healthcare Ltd reserve the right to return product that does not contain a decontamination certificate).
- Mark the package 'Service Department'.

For further details, refer to NHS document HSG(93)26 (UK only).

Service Department,  
Huntleigh Healthcare,  
Diagnostic Products Division,  
35, Portmanmoor Rd.,  
Cardiff. CF24 5HN  
United Kingdom.

Tel: +44 (0)29 20485885  
Fax: +44 (0)29 20492520  
Email: [sales@huntleigh-diagnostics.co.uk](mailto:sales@huntleigh-diagnostics.co.uk)  
[service@huntleigh-diagnostics.co.uk](mailto:service@huntleigh-diagnostics.co.uk)  
Website: [www.huntleigh-diagnostics.com](http://www.huntleigh-diagnostics.com)

### 3.1 Service Life

This is defined as the period during which the device is expected to remain safe and suitable to meet its intended use, and all risk control measures remain effective.

The service life for this device is seven years.

### 3.2 Maintenance and Repair

There are no user serviceable parts inside the Doppler unit or probes. This product does not require periodic maintenance. Inspection is recommended each time the product is used, paying particular attention to the tip of the probes, checking for cracks etc., and to the cable and connector. Any unusual sounds or intermittent behaviour should be investigated.

Spare parts are available. Please refer to service manual for further information and part numbers. A full technical description is provided in the Service Manual 772490.


## Caution

---


**Servicing cannot be carried out while the Doppler is in use.**

## 4. Product Identification


### 4.1 Product Controls

1	Headphone Socket
2	USB Port
3	LCD Panel
4	Function Button 1 / On/Off Button
5	Function Button 2
6	Function Button 3 / Setup
7	Loudspeaker
8	Probe Holder
9	Trolley Mount
10	Volume Up
11	Volume Down
12	Pocket Clip
13	Battery Compartment + Micro SD Card Slot
14	Rear Panel Label


## 4.2 Symbol Identification

	Patient applied parts (ultrasound probes) are type BF according to the definitions in IEC 60601-1.		
	General Warning		Attention, consult accompanying documents / Instructions for Use
	This symbol signifies that this product, including its accessories and consumables is subject to the WEEE (Waste Electrical and Electronic Equipment) regulations and should be disposed of responsibly in accordance with local procedures.		
	This symbol signifies that this product complies with the essential requirements of the Medical Device Directive (93/42/EEC) - Medical Device Regulation (EU/2017/745).		
RxOnly	Federal law restricts this device to sale by, or on the order of a licensed healthcare practitioner.		
Made in the UK by:	<b>Huntleigh Healthcare Ltd.</b> 35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk www.huntleigh-diagnostics.com		
	Legal Manufacturer in association with the CE mark in Europe ArjoHuntleigh AB Hans Michelsensgatan 10 211 20 Malmö, Sweden		
<b>IP20</b>	Protected against ingress of solid foreign objects >12.5 mm diameter. Not protected against ingress of water.		
<b>IPx7</b>	Protected against ingress of water 1 m immersion for up to 30 minutes.		
<b>IPx1</b>	Protected against vertically falling water drops.		
	Power On/Standby		USB Port
	Device Identifier		Serial Number
	Reference Number		Medical Device
	Fragile		Keep Dry
	Atmospheric Pressure Limitations		Relative Humidity Limitations
	Temperature Limitations		Cardboard packaging can be recycled
	LATEX FREE Does not contain Latex		PVC FREE Does not contain PVC
	Headphone Socket		Alignment mark
	Volume Up		Volume Down

**Note: Product labelling should be readable from a distance of up to 0.7m.**

## 4.3 Display Status Bar


The Status Bar appears at the top of most screens, the information displayed depends on the operating mode.

Status Bar Icons			
	Battery Level Low		Date / Time
	USB Connected		USB Transmitting
	USB Receiving		Signal Strength
	Standard Mode - FHR outside user set range		Standard Mode - FHR within user set range
	Smooth Mode - FHR outside user set range		Smooth Mode - FHR within user set range
	Manual Mode		


## 5. Prepare the Doppler for Use

### 5.1 Battery Insertion / Replacement

Disconnect the Doppler from any other equipment before removing the battery cover.

Insert a suitable tool into the recess, release the clip and gently lever off the battery cover.

Insert the batteries according to the diagram, observing polarity.

- Use either alkaline LR6 (non-rechargeable) or NiMH HR6 (rechargeable) batteries.
- Do not mix rechargeable and non-rechargeable batteries.

**Note:** If the Doppler will not be used for an extended period, remove the batteries.

### 5.2 Probe Connection

To connect the probe to the SRX Doppler, align the arrow on the connector with the slot on the probe and push together firmly. (Doppler models SR2 and SR3 have captive probes).




To disconnect the probe, pull the connector from the probe. DO NOT pull the cable.

## 5.3 Change the Doppler Settings

### 5.3.1 Setup Screen


**Note:** A probe must be connected before the Setup Screen can be accessed.

Press  to turn the unit ON, then press and hold button 3 to enter Setup mode.


		Brightness		Date and Time Setup
		Battery type Selection		Lock screen
		FHR Averaging Mode		Trace Options
		Audio Recording		Timer

Press Button 1  to move the selection, Button 2  to accept, or Button 3  to go back.

### 5.3.2 Set the Screen Brightness


### 5.3.3 Set the Date and Time

Press  to accept or  to discard the changes.

### 5.3.4 Battery Type Selection

HR6 = rechargeable NiMH batteries, LR6 = non-rechargeable alkaline batteries

## 6. Operation

Press  and hold for one second to turn the unit ON. The Doppler displays the Obstetric Live Waveform Screen.

### 6.1 SRX Probe options

Two types of probe are available:

OP2XS	2 MHz $\pm$ 1%	OP3XS	3 MHz $\pm$ 1%
-------	----------------	-------	----------------

The Obstetric FHR Numeric Screen is displayed when an obstetric probe is detected. The strength of the Doppler signal is displayed in the status bar at the top of the display.

### 6.2 SR2 & SR3 Waterproof Probes

SR2 is supplied with a captive waterproof 2 MHz probe and cable while SR3 is supplied with a captive waterproof 3 MHz probe and cable.

When used underwater, refer to the Cleaning Instructions in section 7, to prevent cross-contamination.

#### Patient Preparation

Make the patient comfortable in a semi-recumbant or sitting position. Apply a liberal amount of gel to the abdomen. Place the faceplate of the probe flat against the abdomen above the symphysis pubis. Adjust the probe to get an optimum audio signal, ideally by angling the probe around while maintaining firm pressure. Avoid sliding the probe over the skin. Best performance is from the fetal heart itself, characterised by 'slapping' valve sounds rather than umbilical artery or placental sounds. The FHR, averaged over 4-heart beats, is displayed on the 3-digit readout.

**\*Note: Gel is not required when SR2 and SR3 model probes are used underwater.**

### 6.3 Obstetric Monitoring

#### Obstetric FHR Numeric Mode Screen

In this Mode, the FHR is displayed (in beats per minute) in large digits and is continually updated. When the rate cannot be determined, 3 x dashes are displayed.

Press  to switch to FHR Trace Display.

### 6.3.1 Obstetric FHR Trace Mode Screen\*

Trace display mode is for indication only and is not a replacement for conventional fetal monitoring. If any concern arises from viewing this trace, alternative means (e.g. full CTG) must be used to determine fetal condition.

This screen displays the FHR as a trace on a graph. The horizontal and vertical scales are determined in setup\*\*. The trace displays from left to right until it fills the screen then scrolls to the left as each new data point is added.

FHR Trace Mode Icons	
	Numeric Display Mode
	Register fetal movement and place marker
	Stop Trace and display Frozen Trace option

\*\* Trace scaling is equivalent to 1 or 3 cm/min and 20 or 30 beats per cm, scaled down to screen size, maintaining the same aspect ratio to avoid distortion of the trace for easy visual interpretation.

\*Note: These options are only available if the FHR Trace Mode screen has been enabled.

### 6.3.2 Trace Scroll section

Press  when the desired information is displayed.


## 6.4 Frozen Trace

Press  at any time to return to the live trace screen.

To Save a trace:

To open a saved trace:

## 6.5 Battery Charging

- Only HR6 (NiMH) batteries can be charged. Check the battery type before connecting the charger.
- Do not attempt to recharge normal alkaline batteries. They may leak, cause a fire or even explode.
- Only use the charger and USB lead supplied by Huntleigh.
- Switch off the Doppler before charging.
- Do not use the Doppler on patients when the charger is connected.

- Insert the supplied charger lead into the USB socket at the top of the Doppler.
- Connect the charger and switch on mains power.

Charging will take approximately 5 - 6 hours, depending on battery status.

If the battery level is low, a  symbol will appear on the Status bar.

During charging, a  symbol will be displayed on the screen.

When fully charged, the symbol will change to .

## 6.6 Data transfer to a PC

The Doppler does not contain patient information however, stored traces and data can be transferred to a PC running Huntleigh reporting software, via the same USB socket as for battery charging. Consult your local sales representative for details.

**Note:** Please observe the safety warnings in the Huntleigh software IFU.

## 6.7 After Use

Press and hold the On/Off button for one second to switch the unit off.

Refer to the cleaning section before storing or using the unit on another patient.

# 7. Care and Cleaning

## 7.1 General Care

The Doppler contains delicate components, for example the probe tip, which should be handled and treated with care. Periodically, and whenever the integrity of the system is in doubt, carry out a check of all functions as described in the relevant section of this IFU. If there are any defects, contact Huntleigh or your distributor for repair or to order a replacement.

## Cautions

- Check with your facility's local infection control policy and medical equipment cleaning procedures.
- Observe warnings and guidance on cleaning fluid labelling regarding use and personal protective equipment (PPE).
- If detergent or disinfectant wipes are used ensure that excess solution is squeezed from the wipe prior to use.
- Always switch off the Doppler and disconnect from the AC supply before cleaning and disinfecting.
- Always wipe off disinfectant using a cloth dampened with clean water.
- Do not allow any fluid to enter the products and do not immerse in any solution.
- Do not use abrasive cloths or cleaners.
- Do not use automatic washers or autoclaves.
- Do not use Phenolic detergent based disinfectants, solutions containing cationic surfactants, ammonia based compounds or perfumes and antiseptic solutions.

## 7.2 Cleaning and Disinfecting the Doppler

Always keep the external surfaces clean and free of dirt and fluids using a clean dry cloth.

- **Wipe any fluids from the surface of the product using a clean dry cloth.**
- **Wipe with a cloth dampened in 70% Isopropyl Alcohol.**
- **Completely dry with a clean, dry cloth.**
- **If the product has been contaminated use the methods described for probes.**

## 7.3 Cleaning and Disinfecting Probes

Clean the probes before examining a patient using low risk cleaning method below.

Following patient examination, clean and/or disinfect the probes by the appropriate method based upon the level of cross contamination risk, as defined below:

Risk	Definitions	Procedure
Low	Normal use or low risk situations include patients having intact skin and no known infection.	1. Remove soiling, wipe with a mild neutral detergent and then wipe with a cloth dampened in water. 2. Completely dry with a clean cloth.
Medium	The patient has a known infection, skin is not intact, the part is heavily soiled.	1. Follow low risk procedure then wipe with a cloth dampened in Sodium Hypochlorite (1,000ppm). 2. After two minutes wipe with a cloth dampened in water and then dry with a clean cloth.
High	This procedure should only be used when the part has been contaminated by blood.	1. Follow low risk procedure then wipe with a cloth dampened in Sodium Hypochlorite (10,000ppm). 2. After two minutes wipe with a cloth dampened in water and then dry with a clean cloth.

### Caution

**Repeated and unnecessary use of concentrated solutions will result in damage to the product. Do not allow Sodium Hypochlorite solutions to come into contact with metal parts.**

The use of disinfectant materials other than those listed is the responsibility of the user for their efficacy and compatibility with the device.

## 8. Troubleshooting

This section lists some of the more common problems encountered during use together with possible causes. If the problem cannot be located after consulting the table in this section, turn off the Doppler and consult a qualified technician. Before attempting trouble-shooting, verify that the batteries are charged.

SYMPTOM	POSSIBLE CAUSE / REMEDY
Doppler will not turn on.	Replace / recharge batteries.
Audio Only	Doppler model does not support visual functionality
No Audio signal	Incorrect volume setting
Poor Signal	Probe / sensor incorrectly positioned, or Insufficient Gel
No Signal	Damaged probe / sensor, or Incorrect probe / sensor
Screen displays:	 Damaged probe / sensor, or No Probe
	 Incompatible probe / sensor, or Incorrect Probe / sensor
	 Incorrect battery fitted

## 9. Specifications

### 9.1 Equipment Classification


Type of protection against electric shock.	Internally powered equipment
Degree of protection against electric shock 	Type BF - equipment with an applied part 
Mode of operation.	Continuous
Degree of protection against harmful ingress of particles and/or water.	Main Unit: IP20*, SR2/SR3 probes : IPX7 Protective Pouch : IP22
Degree of safety of application in the presence of a flammable anaesthetic	Equipment not suitable for use in the presence of a FLAMMABLE ANAESTHETIC MIXTURE WITH AIR, OXYGEN OR NITROUS OXIDE

\*For home use, this can be upgraded to IPx2 when using the protective pouch (ACC-OBS-080).

### 9.2 FHR Performance

Standard Mode	Range: 60 - 210 bpm Averaging: 4 beats	Resolution: 1 bpm Accuracy: $\pm 3$ bpm
Smoothed Mode	Range: 60 - 210 bpm Averaging: 8 beats	Resolution: 1 bpm Accuracy: $\pm 3$ bpm
Manual Mode	Range: 60 - 210 bpm Averaging: 10 beats	Resolution: 1 bpm Accuracy: $\pm 3$ bpm

### 9.3 General

Charger- 'R' models only (Part No. 772559)	Protection : Input Voltage : Output Voltage : Input Frequency : Standby power consumption:	Class II  100-240 V AC $\pm 10\%$ 5VDC $\pm 5\%$ 50 - 60Hz 230V AC $\leq 0.1W$
Max. Audio Output	500 mW rms typical (loudspeaker)	
Auto shut-off	3 minutes	
Headphone output	Max. output Power: Connector:	25 mW rms (32 $\Omega$ ) 3.5 mm stereo jack socket
USB Port	Micro USB	<b>SD Card Slot</b> Micro SD
Real Time clock battery	RENATA CR1025, 3V Lithium	
Battery Type	LR6 (Alkaline cells 1.5V) or HR6 (NiMH rechargeable 1.2V)	
Battery Life	Typically, 500 x 1 minute examinations Note: Battery life is typically 2 years or 500 charge/discharge cycles	
Software Version	v4.0	
Size	140 x 33 x 75 mm	<b>Weight</b> 280 g

## 9.4 Environmental

Operating	
Temperature range	+5°C to +40°C
Relative Humidity	15% to 90% (non condensing)
Pressure	700 hPa to 1060 hPa
Transport and Storage between uses	
Without relative humidity control	-25°C to +5°C
At a r.h. of up to 90% non-condensing	+5°C to +35°C
At a water vapour pressure up to 50 hPa	>+35°C to +70°C

## 9.5 Standards Compliance

IEC 60601-1: 2012	IEC 60601-1-11: 2015
EN 60601-2-37: 2015 Thermal Indices (TI) and Mechanical Index (MI) are below 1.0 for all device settings.	IEC 60601-1-2: 2014

## 9.6 Accessories

Use only the recommended accessories. See [www.huntleigh-diagnostics.com](http://www.huntleigh-diagnostics.com) for a list of accessories.

# 10. Electromagnetic Compatibility

Make sure the environment in which the Doppler is installed is not subject to strong sources of electromagnetic interference (e.g. radio transmitters, mobile phones). This equipment generates and uses radio frequency energy. If not installed and used properly, in strict accordance with the manufacturer's instructions, it may cause or be subject to interference. Type-tested in a fully configured system, complies with IEC 60601-1-2, the standard intended to provide reasonable protection against such interference. Whether the equipment causes interference may be determined by turning the equipment off and on. If it does cause or is affected by interference, one or more of the following measures may correct the interference:


- Reorient the equipment
- Relocate the equipment with respect to the source of interference
- Move the equipment away from the device with which it is interfering

## Warnings

- **The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of the Doppler as replacement parts for internal components, may result in increased emissions or decreased immunity of the Doppler.**
- **The Doppler should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the Doppler should be observed to verify normal operation in the configuration in which it will be used.**
- **Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Doppler including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.**

This section is only applicable to United Kingdom (UK) market when UK marking is applied to the Arjo medical device labelling.

**UK Symbol:**

**UK Responsible Person & UK Importer:**

**Arjo (UK) Ltd., ArjoHuntleigh House, Houghton Regis. LU5 5XF**

Is the appointed UK Responsible Person as defined in UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended).

For Northern Ireland (NI) CE marking will still apply until further amendment to applicable regulations.

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This page is intentionally left blank.

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Manufactured in the UK by **Huntleigh Healthcare Ltd** on behalf of;

ArjoHuntleigh AB  
Hans Michelsensgatan 10  
211 20 Malmö, Sweden

2797

**Huntleigh Healthcare Ltd.**

35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom  
T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk  
www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office:  
ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF  
©Huntleigh Healthcare Limited 2021

**A Member of the Arjo Family**

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

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[www.huntleigh-diagnostics.com/](http://www.huntleigh-diagnostics.com/)

[www.huntleigh-healthcare.us/](http://www.huntleigh-healthcare.us/)

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**Huntleigh Healthcare GmbH**