

# HUNTLEIGH

A MEMBER OF THE ARJO FAMILY



## Dopplex

Vascular Assessment  
Accessories & Consumables  
Catalogue

2023



# CONTENTS

## **DOPPLER ACCESSORIES** **A**

High Sensitivity Probes	<b>3</b>
Intraoperative Probes / Packs	<b>4</b>
Generic Accessories	<b>5</b>
Generic Consumables	<b>9</b>

## **ABILITY ACCESSORIES** **B**

Cuff Kits	<b>10</b>
Consumables	<b>12</b>

## **LEGACY ACCESSORIES** **C**

Dopplers	<b>13</b>
----------	-----------



**VP4XS**

High Sensitivity 4MHz probe  
(Detection of Deep Lying vessels)



**VP5XS**

High Sensitivity 5MHz probe  
(Oedematous Limbs & Deep Lying Vessels)



**VP8XS**

High Sensitivity 8MHz probe  
(Detection of Peripheral Vessels & Calcified Arteries)



**VP10XS**

High Sensitivity 10MHz probe  
(Detection of Smaller Vessels)



**EZ8XS**

High Sensitivity Wide-beam 8MHz probe  
(Detection of Peripheral Vessels & Calcified Arteries)



**PPGA1**

Arterial PPG Adaptor and Probe  
(Detection of Toe Pressures - Used in ATP Kits)



 Compatible

 Non-Compatible

DMX  
SD2  
D900

**DIPP10**

Pre-sterilised Single-Use Intraoperative 8MHz Doppler probes (Box of 10)



**DISP10XS**

Single-Use Intraoperative Probe Starter Pack: Includes PA8XS Probe Adaptor, Pole Clamp & DIPP10



**DISP10-DMXR**

Single-Use Intraoperative Probe Starter Pack: Includes PA8XS Probe Adaptor, Pole Clamp & Box of 10 Single-use Intraoperative 8MHz Doppler Probes and DMXR Doppler



**PA8XS**

Intraoperative Probe Adapter - DMX range



**ISP-XS**

Intraoperative Probe Starter Pack: Includes PA8XS Probe Adaptor, Clamp



**ACC47**

IV Pole Clamp for use with Intraoperative Probe



**ACC21**

Stereo Headset



**ACC34**

Soft Carry Pouch



**ACC232**

Dopplex Stand - shipped flat packed, Includes Adaptor Plate



**ACC159**

Adult Arm/Ankle Cuff with Colder connector (23-33cm)



**ACC160**

Large Arm/Ankle Cuff with Colder connector (31-40cm)



**ACC161**

Thigh Cuff with Colder Connector (38-50cm)



Compatible

Non-Compatible

**ACC162**

Large Digit/Toe Cuff with Colder Connector  
(2.5cm) (pack of 5)

**ACC163**

Standard Digit/Toe Cuff with Colder Connector  
(1.9cm) (pack of 5)

**ACC185**

Digit /Toe Cuff with colder connector  
- 3 Large & 3 Standard

**ACC164**

Tourniquet Cuff with Colder Connector

**ACC168**

Aneroid Sphygmomanometer with trigger  
release valve and Colder Connector

**ACC187**

Tourniquet cuff & Aneroid Sphygmomanometer  
Kit with trigger release valve & Colder Connector

**ACC188**

Blood Pressure cuff & Aneroid  
Sphygmomanometer Kit with trigger release  
valve & Colder Connector

DMX  
SD2  
D900



**ACC226**

Medical Grade Rechargeable Kit



**ACC227**

Micro SD card



**ACC228**

Doppler Stand Adaptor plate



**ACC-VAS-028**

T-Tube Toe Cuff Inflator System



**ACC-VAS-029**

Large Carry bag - use with ATP/ABI Kits



Compatible

Non-Compatible



DMX  
SD2  
D900



**ACC-VAS-030**

USB cable use with Vascular Reporter (3m)



**ACC-OBS-080**

Splash-proof Doppler pouch (IPX2)



**ACC3**

Ultrasound Gel  
(Box 12 x 250ml Bottles)



**ACC24**

Ultrasound Gel  
(Box 12 x 60ml Tubes)



**ACC-1300-0145**

Ultrasound Gel (box 48 x 20gm Sachet)





**ACC179**  
APPG Adhesive Straps  
(Pack of 240)



**ACC224**  
Neuropen & Tip  
- Includes 10g Monofilament



# DOPPLEX ABILITY



High Sensitivity Probes

Intraoperative Probes / Packs

Generic Accessories

Generic Consumables

## ACC-VAS-007

Large Adult Right Arm Cuff - Red tube  
(34-46cm)

**RIGHT**  
**LEFT**  
**KIT**



## ACC-VAS-008

Large Adult Right Ankle Cuff - Black tube  
(24-35cm)



## ACC-VAS-009

Large Adult Left Arm Cuff - Yellow tube  
(34-46cm)



## ACC-VAS-010

Large Adult Left Ankle Cuff - Green tube  
(24-35cm)



**B**

■ Compatible

□ Non-Compatible

**ACC-VAS-011**

Large Adult Cuff Set  
(34-46cm)

*RIGHT*  
*LEFT*  
*KIT*



**ACC-VAS-023**

Adult Right Arm Cuff - Red tube  
(22-36cm)



**ACC-VAS-024**

Adult Right Ankle Cuff - Black tube  
(18-28cm)



**ACC-VAS-025**

Adult Left Arm Cuff - Yellow tube  
(22-36cm)



**ACC-VAS-026**

Adult Left Ankle Cuff - Green tube  
(18-28cm)



**ACC-VAS-027**

Adult Arm Cuff Set  
(22-36cm)



Compatible

Non-Compatible

**ACC-VAS-012**

Fixing Plate (required for ACC-VAS-013  
or ACC-VAS-154)

**ACC-VAS-013**

Pole Stand / Trolley  
(Requires ACC-VAS-012)

**ACC-VAS-015**

Carry bag

**ACC-VAS-154**

Wall Mount (requires ACC-VAS-012)

**ACC-VAS-016**

Infection Control Sleeves  
(Box of 100)

**ACC-VAS-017**

Standard Thermal Paper - Pack of 5 (approx.  
55 prints per roll)

**ACC-VAS-019**

Self-adhesive Thermal Paper - Pack of 5  
(approx. 55 prints per roll)

# LEGACY



## ACC171

Rechargeable Battery Kit (UK)

## ACC172

Rechargeable Battery Kit (USA)

## ACC173

Rechargeable Battery Kit (EU)

-Supplied with Charger & 2 Batteries

## PA8-HG

Intraoperative Probe Adaptor

## DISP10

Single-Use Intraoperative Probe Starter Pack:  
Includes PA8-HG Probe Adaptor, Pole Clamp  
& Box of 10 Single-use Intraoperative 8MHz  
Doppler Probes

## ISP

Intraoperative Probe Starter Pack:  
Includes PA8-HG Probe Adaptor, Clamp

D900  
SD2  
MD200



Compatible

Non-Compatible

**ACC232**

Dopplex Stand - Shipped flat packed



**ACC180**

Dopplex Stand - Shipped flat packed



**ACC190**

RS232-USB Adaptor for use with Dopplex Reporter Software Package (DR4), running on a PC without a RS232 port - requires ACC35 Cable





As a proud member of the Arjo family, we have been committed to supporting healthcare professionals in improving outcomes and enhancing patient wellbeing since 1979. We do this through our proven solutions for Vascular Assessment & Treatment and Fetal & Patient Monitoring. With innovation and customer satisfaction as our guiding principles, we strive for clinical excellence and improved performance, for life.

**Manufactured and distributed by Huntleigh Healthcare Ltd.**  
35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom  
T: +44 (0)29 20485885 | sales@huntleigh-diagnostics.co.uk  
www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office:  
Arjo Huntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF  
© Huntleigh Healthcare Limited 2019

**A Member of the Arjo Family**

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

4401000024

**HUNTLEIGH**

LIT 710507-12